


# Matching accessoires: Immersion sleeves

Brass immersion sleeve R1/2 " TH08-MS				
Article no.	809521 00XX	Code		ML (mm)
Temp. max	+150 °C	05		50
Pressure proof up to	10 bar (statisch)	10		100
Quick-release strap	Nickel plated	15		150
Process connection	R1/2 "	20		200
Wrench size	22	25		250
Ø Inside / Outside (mm)	6,5 / 8	30		300
Scope of delivery	Brass immersion sleeve R1/2 " TH-MS	40		400
Order code	809521 00	--		



Please attach the code for the mounting length.

### Technical drawing




Stainless steel immersion sleeve G1/2 " TH08-VA				
Article no.	809521 10XX	Code		ML (mm)
Temp. max	+600 °C	05		50
Pressure proof up to	40 bar (statisch)	10		100
Quick-release strap	1.4571 (316TI)	15		150
Process connection	G1/2 "	20		200
Wrench size	27	25		250
Ø Inside / Outside (mm)	6,5 / 8	30		300
Scope of delivery	Stainless steel immersion sleeve G1/2 " TH-VA	40		400
Your order code	809521 10	--		



Please attach the code for the mounting length.

### Technical drawing



MW / KS / 22.11.2023

### Testo Sensor GmbH

Testo-Straße 1  
D-79853 Lenzkirch

+49 7653 96597-0  
+49 7653 96597-99

info@testo-sensor.de  
www.testo-sensor.de

You can find our standard portfolio in our  
webshop at: [www.testo-sensor.shop](http://www.testo-sensor.shop)


**Stainless steel immersion sleeve G1/2 " with neck tube TH08-VA**



Article no.	809521 20XX	Code	EL (mm)
Temp. max	+600 °C	05	50
Pressure proof up to	40 bar (statisch)	10	100
Quick-release strap	1.4571 (316TI)	15	150
Process connection	G1/2 "	20	200
Wrench size	27	25	250
Ø Inside / Outside (mm)	6,5 / 8	30	300
Length neck tube (mm)	90		
Scope of delivery	Stainless steel immersion sleeve G1/2 " with neck tube TH-VA		
Your order code	809521 20	--	

Please attach the code for the mounting length.

**Technical drawing**



MW / KS / 22.11.2023

**Testo Sensor GmbH**

Testo-Straße 1  
D-79853 Lenzkirch

+49 7653 96597-0  
+49 7653 96597-99

info@testo-sensor.de  
www.testo-sensor.de

You can find our standard portfolio in our  
webshop at: [www.testo-sensor.shop](http://www.testo-sensor.shop)